


Mirror Soft® Atraumatic Blunt Cannula

Atraumatic Blunt Cannula by Mirror Soft® for
Biostimulators and Dermal Fillers



Less Cannula Changes, More Safety and Patient Comfort

An extra wide exit opening and polished inner tube ensuring unhindered flow of suspension based biostimulators and dermal fillers and reducing the risk of clogging and cannula change during injection procedures.



Safe and Practical

The pilot needle and atraumatic blunt cannula come in one sterile blister pack. Less additional handling steps reduce chances of infection and enhance patient safety.



Innovation That Pays

The efficiency in our state-of-the-art manufacturing facilities allows us to cut costs dramatically while maintaining the highest standards set by the European EMA. We are proud to be able to forward the cost efficiency of our innovations to our customers offering Mirror Soft® Blunt Cannulas at a price significantly below market average.

Size	Outer Diameter	Length	Length	Length
17G	1,40mm	-	50mm	70mm
18G	1,20mm	-	50mm	70mm
19G	1,07mm	-	50mm	70mm
20G	0,90mm	-	50mm	70mm
21G	0,80mm	-	50mm	70mm
22G	0,70mm	-	50mm	70mm
23G	0,60mm	30mm	50mm	70mm
24G	0,56mm	-	50mm	-
25G	0,50mm	40mm	50mm	70mm
27G	0,41mm	30mm	40mm	50mm
30G	0,31mm	25mm	-	-

Visit us and order here:



The information provided in this data sheet is for general guidance only and may be subject to change. It is the user's responsibility to ensure the product is suitable for its intended purpose and used according to the manufacturer's instructions. Single use only. Do not reuse or re-sterilize.

Promainject is part of PromaMedical GmbH
Wernher-von-Braun-Straße 10
85640 Putzbrunn – Munich,
Germany

T: +49 (0)89 7474713-0
F: +49 (0)89 7474713-1
M: info@promainject.com
W: www.promainject.com

e-mail us:

