



Proma inject LDS Syringe 1mL



Low-Dead-Space (LDS) Luer-Lock 1mL Syringe by Proma inject for Aesthetic Applications



More Product for the Patient


The Low-Dead-Space technology based on an exact fit between rubber stopper and needle hub allows Prima inject LDS Syringes close to zero product loss. The efficiency is even higher when combined with a low-dead-space needle.



Safety and Comfort

The luer-lock adapter tightly secures the injection needle on the syringe allowing the injector to work with a wider range of pressure on the plunger without risking the needle to detach.

The specially designed plunger slides without blocking and ensures high level of control during the injection.



Innovation That Pays

The efficiency in our state-of-the-art manufacturing facilities allows us to cut costs dramatically while maintaining the highest standards set by the American FDA and European EMA.

We are proud to be able to forward the cost efficiency of our innovations to our customers offering the Proma inject LDS Syringes at a price significantly below market average.

Visit us and order here:



The information provided in this data sheet is for general guidance only and may be subject to change. It is the user's responsibility to ensure the product is suitable for its intended purpose and used according to the manufacturer's instructions. Single use only. Do not reuse or re-sterilize.

PromalInject is part of PromaMedical GmbH
Wernher-von-Braun-Straße 10
85640 Putzbrunn – Munich,
Germany

T: +49 (0)89 7474713-0
F: +49 (0)89 7474713-1
M: info@promainject.com
W: www.promainject.com

e-mail us:

