


ZELOSTAT™ Ultra-Low-Dead-Space Injection Needle




Ultra-Low-Dead-Space (ULDS) Wide-Diameter Hypodermic Needle by ZELOSTAT™ for Botulinum Toxin Micro-Injections



More Product for Your Patient

The Low-Dead-Space technology based on a soft polymer ring filling any extra space, dramatically reducing product waste.

Changing the needles 3 times will already save over 100µL of product or one extra injection point you can provide to your patient. The efficiency is even higher when combined with a Low-Dead-Space Syringe by Proma inject.




Keep Your Routine

With an inner diameter of 0.1mm, our needle offers a flow and plunger pressure comparable to a conventional 34G, allowing you to maintain your current injection technique. The injector feels the same, but the patient experiences less needle.

More Patient Comfort, Less Bruising

Approx. 25% thinner than conventional 34G needles. Our 35G ULDS Needle reduces the risk of adverse events and elevates patient comfort during the injection procedure – less visible than invisible, it's ethereal.



Article Number	Outer Diameter	Length	Packaging
100533	31G (0.25mm)	13mm	100pcs
100532	31G (0.25mm)	9mm	100pcs
100525	31G (0.25mm)	4mm	100pcs
100539	34G (0.18mm)	9mm	100pcs
100538	34G (0.18mm)	4mm	100pcs
100543	35G (0.16mm)	4mm	100pcs
100543-1	35G (0.16mm)	4mm	50pcs

Visit us and order here:



The information provided in this data sheet is for general guidance only and may be subject to change. It is the user's responsibility to ensure the product is suitable for its intended purpose and used according to the manufacturer's instructions. Single use only. Do not reuse or re-sterilize.

Promainject is part of PromaMedical GmbH
Wernher-von-Braun-Straße 10
85640 Putzbrunn – Munich,
Germany

T: +49 (0)89 7474713-0
F: +49 (0)89 7474713-1
M: info@promainject.com
W: www.promainject.com

e-mail us:

